

Product Data Sheet

Recombinant Human DAPP1 Protein

Catalog #	Source	Reactivity	Applications
CST1833	E. coli	H	
Description	Recombinant Human DAPP1 Protein is produced by E. coli expression system and the target gene encoding Met11-Thr160 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	DAPP1		
Alternative Names	BAM32; Dual adapter for phosphotyrosine and 3-phosphotyrosine and 3-phosphoinositide; hDAPP1; B lymphocyte adapter protein Bam32; B-cell adapter molecule of 32 kDa		
Entrez Gene	27071		
SwissProt	Q9UN19		
Purity			
Chemical Structure	MSTQDPSDLW SRSDGEAELL QDLGWYHGNI TRHAAEALL SNGCDGSYLL RDSNETTGLY SLSVRAKDSV KHFFHVEYTG YSFKFGFNEFS SLKDFVKHFA NQPLIGSETG TLMVLKHPYP RKVEEPSIYE SVRVHTAMQT GRTEDDLVP		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com