

Product Data Sheet

Recombinant Human KAZALD1 Protein

Catalog #	Source	Reactivity	Applications
CST1847	E. coli	H	
Description	Recombinant Human KAZALD1 Protein is produced by E. coli expression system and the target gene encoding Thr181-Asn300 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	KAZALD1		
Alternative Names	Kazal-type serine protease inhibitor domain-containing protein 1		
Entrez Gene	81621		
SwissProt	Q96I82		
Purity			
Chemical Structure	TWNVTGQDVI FGCEVFAYPM ASIEWRKDGL DIQLPGDDPH ISVQFRGGPQ RFEVTGWLQI QAVRPSDEGT YRCLGRNALG QVEAPASLTV LTPDQLNSTG IPQLRSLNLV PEEEAEESEN		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com