

Product Data Sheet

Recombinant Human RNF31 Protein

Catalog #	Source	Reactivity	Applications
CST1849	E. coli	H	
Description	Recombinant Human RNF31 Protein is produced by E. coli expression system and the target gene encoding Val971-Arg1070 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	RNF31		
Alternative Names	ZIBRA; E3 ubiquitin-protein ligase RNF31; HOIL-1-interacting protein; HOIP; RING finger protein 31; Zinc in-between-RING-finger ubiquitin-associated domain protein		
Entrez Gene	55072		
SwissProt	Q96EP0		
Purity			
Chemical Structure	VIEQKEVPNG LRDEACGKET PAGYAGLCQA HYKEYLVSLI NAHSLDPATL YEVEELETAT ERYLHVPRQP LAGEDPPAYQ ARLLQKLTEE VPLGQSIPRR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com