

Recombinant Human FAM3C Protein

Catalog #	Source	Reactivity	Applications
CST1851	E. coli	H	
Description	Recombinant Human FAM3C Protein is produced by E. coli expression system and the target gene encoding Gln25-Asp227 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FAM3C		
Alternative Names	ILEI; Protein FAM3C; Interleukin-like EMT inducer		
Entrez Gene	10447		
SwissProt	Q92520		
Purity			
Chemical Structure	QVFEIKMDAS LGNLFARSAL DTAARSTKPP RYKCGISKAC PEKHFAFKMA SGAANVVGPK ICLEDNVLMG GVKNNVGRGI NVALANGKTG EVLDTKYFDM WGGDVAPFIE FLKAIQDGTI VLMGTYDDGA TKLNDEARRL IADLGSTSIT NLGFRDNWVF CGGKGIKTKS PFEQHIKNNK DTNKYEGWPE VVEMEGCIPQ KQD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com