

Product Data Sheet

Recombinant Human MCFD2 Protein

Catalog #	Source	Reactivity	Applications
CST1856	E. coli	H	
Description	Recombinant Human MCFD2 Protein is produced by E. coli expression system and the target gene encoding Glu27-Gln146 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	MCFD2		
Alternative Names	SDNSF; Multiple coagulation factor deficiency protein 2; Neural stem cell-derived neuronal survival protein		
Entrez Gene	90411		
SwissProt	Q8NI22		
Purity			
Chemical Structure	EEPAASFSQP GSMGLDKNTV HDQEHIMEHL EGVINKPEAE MSPQELQLHY FKMHDYDGNN LLDGLELSTA ITHVHKEEGS EQAPLMSEDE LINIIDGVLR DDDKNNDGYI DYAEFAKSLQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com