

Recombinant Human UNC5D Protein

Catalog #	Source	Reactivity	Applications
CST1862	E. coli	H	
Description	Recombinant Human UNC5D Protein is produced by E. coli expression system and the target gene encoding Ala41-Arg240 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	UNC5D		
Alternative Names	KIAA1777; UNC5H4; Netrin receptor UNC5D; Protein unc-5 homolog 4; Protein unc-5 homolog D		
Entrez Gene	137970		
SwissProt	Q6UXZ4		
Purity			
Chemical Structure	ALPESIPSAP GTLPHFIEEP DDAYIIKSNP IALRCKARPA MQIFFKCNGE WVHQNEHVSE ETLDESSGLK VREVFINVTR QQVEDFHGPE DYWCQCVAWS HLGTSKSRKA SVRIAYLRKN FEQDPQGREV PIEGMIVLHC RPPEGVPAAE VEWLKNEEPI DSEQDENIDT RADHNLIRQ ARLSDSGNYT CMAANIVAKR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com