

Product Data Sheet

Recombinant Human FABP9 Protein

Catalog #	Source	Reactivity	Applications
CST1875	E. coli	H	
Description	Recombinant Human FABP9 Protein is produced by E. coli expression system and the target gene encoding Met1-Val132 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FABP9		
Alternative Names	Fatty acid-binding protein 9; Testis lipid-binding protein; TLBP; Testis-type fatty acid-binding protein; T-FABP		
Entrez Gene	646480		
SwissProt	Q0Z7S8		
Purity			
Chemical Structure	MVEPFLGTWK LVSENFEDY MKELGVNFAA RNMAGLVKPT VTISVDGKMM TIRTESSFQD TKISFKLGEE FDETTADNRK VKSTITLENG SMIHVQKWLG KETTIKRKIV DEKMVVECKM NNIVSTRIYE KV		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com