

Product Data Sheet

Recombinant Human IGSF3 Protein

Catalog #	Source	Reactivity	Applications
CST1936	E. coli	H	
Description	Recombinant Human IGSF3 Protein is produced by E. coli expression system and the target gene encoding Ala951-Glu1100 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	IGSF3		
Alternative Names	EWI3; KIAA0466; Immunoglobulin superfamily member 3; IgSF3; Glu-Trp-Ile EWI motif-containing protein 3; EWI-3		
Entrez Gene	3321		
SwissProt	O75054		
Purity			
Chemical Structure	ASLQVDTVVP NATVSEKAAF QLDCSIVSRS SQDSRFAVAW YSLRTKAGGK RSSPGLEEQE EEREEEEED DDDDDDDPTER TALLSVGPDA VFGPEGSPWE GRLRFQRLSP VLYRLTVLQA SPQDTGNYSC HVEEWLPSPQ KEWYRLTEEE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com