

# Product Data Sheet

## Recombinant Human CHAD Protein

| Catalog #                 | Source                                                                                                                                                                     | Reactivity | Applications |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|--------------|
| CST1939                   | E. coli                                                                                                                                                                    | H          |              |
| <b>Description</b>        | Recombinant Human CHAD Protein is produced by E. coli expression system and the target gene encoding Asp261-Thr350 is expressed with N-terminus SUMO & C-terminus His tag. |            |              |
| <b>Form</b>               | Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.                                                                                                                       |            |              |
| <b>Gene Symbol</b>        | CHAD                                                                                                                                                                       |            |              |
| <b>Alternative Names</b>  | SLRR4A; Chondroadherin; Cartilage leucine-rich protein                                                                                                                     |            |              |
| <b>Entrez Gene</b>        | 1101                                                                                                                                                                       |            |              |
| <b>SwissProt</b>          | O15335                                                                                                                                                                     |            |              |
| <b>Purity</b>             |                                                                                                                                                                            |            |              |
| <b>Chemical Structure</b> | DGAFLGVTTL KHVHLENNRL NQLPSNFPFD SLETLALTNN PWKCTCQLRG LRRWLEAKAS<br>RPDATCASPA KFKGQHIRT DAFRSCKFPT                                                                       |            |              |
| <b>Quality Control</b>    |                                                                                                                                                                            |            |              |
| <b>Directions for Use</b> | Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.                                                          |            |              |
| <b>Storage/Stability</b>  | Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.                           |            |              |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com