

Product Data Sheet

Recombinant Human ITGB1BP1 Protein

Catalog #	Source	Reactivity	Applications
CST1941	E. coli	H	
Description	Recombinant Human ITGB1BP1 Protein is produced by E. coli expression system and the target gene encoding Thr21-Pro200 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ITGB1BP1		
Alternative Names	ICAP1; Integrin beta-1-binding protein 1; Integrin cytoplasmic domain-associated protein 1; ICAP-1		
Entrez Gene	9270		
SwissProt	O14713		
Purity			
Chemical Structure	TKSKVDSSL GGLSRSTVA SLDTDSTKSS GQSNNSDTC AEFRIKYVGA IEKCLKSEGG GLEGPLDLIN YIDVAQQDGK LPFVPPEEEF IMGVSKYGIK VSTSDQYDVL HRHALYLIIR MVCYDDGLGA GKSLALKTT DASNEEYSLW VYQCNSLEQA QAICKVLSTA FDSVLTSEKP		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com