

Product Data Sheet

Recombinant Human AMN Protein

Catalog #	Source	Reactivity	Applications
CST1947	E. coli	H	
Description	Recombinant Human AMN Protein is produced by E. coli expression system and the target gene encoding Asp131-Gly340 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	AMN		
Alternative Names	Protein amnionless		
Entrez Gene	81693		
SwissProt	Q9BXJ7		
Purity			
Chemical Structure	DAERVPCRHD DVFFPPSASF RVGLGPGASP VRVRSISALG RTFTRDEDLA VFLASRAGRL RFHGPGALSV GPEDCADPSG CVCGNAAEQP WICAALLQPL GGRCPQAACH SALRPQGQCC DLCGAVVLLT HGPAFDLERY RARILDFTLG LPQYHGLQVA VSKVPRSSRL READTEIQVV LVENGPETGG AGRLARALLA DVAENGEALG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com