

Product Data Sheet

Recombinant Human BCL2L11 Protein

Catalog #	Source	Reactivity	Applications
CST1950	E. coli	H	
Description	Recombinant Human BCL2L11 Protein is produced by E. coli expression system and the target gene encoding Glu11-Arg120 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	BCL2L11		
Alternative Names	BIM; Bcl-2-like protein 11; Bcl2-L-11; Bcl2-interacting mediator of cell death		
Entrez Gene	10018		
SwissProt	O43521-2		
Purity			
Chemical Structure	ECDREGRLQ PAERPPQLRP GAPTSLQTEP QDRSPAPMSC DKSTQTPSP CQAFNHLYSA MASMRQAEP DMRPEIWIAQ ELRRIGDEFN AYYARRVFLN NYQAAEDHPR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com