

Product Data Sheet

Recombinant Human LRP2BP Protein

Catalog #	Source	Reactivity	Applications
CST1958	E. coli	H	
Description	Recombinant Human LRP2BP Protein is produced by E. coli expression system and the target gene encoding Asn41-Ala110 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	LRP2BP		
Alternative Names	KIAA1325; LRP2-binding protein; Megalin-binding protein; MegBP		
Entrez Gene	55805		
SwissProt	Q9P2M1		
Purity			
Chemical Structure	NLVDKALQLL KERILKGDTL AYFLRGQLYF EEGWYEEALE QFEEIKEKDH QATYQLGVMY YDGLGTTLDA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com