

Product Data Sheet

Recombinant Human CHRDL1 Protein

Catalog #	Source	Reactivity	Applications
CST1969	HEK293	H	
Description	Recombinant Human CHRDL1 Protein is produced by HEK293 expression system and the target gene encoding Val261-Arg450 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CHRDL1		
Alternative Names	NRLN1; Chordin-like protein 1; Neuralin-1; Neurogenesin-1; Ventroptin		
Entrez Gene	91851		
SwissProt	Q9BU40		
Purity			
Chemical Structure	VSNGKTYSHG ESWHPNLRAF GIVECVLCTC NVTKQECKKI HCPNRYPCY PQQIDGKCKK VCPGKKAKEL PGQSFDNKG YFCGEETMPVY ESVMEDGET TRKIALETER PPQVEVHVWT IRKGILQHFH IEKISKRMFE ELPHFKLVTR TTLSQWKIFT EGQAQISQMC SSRVCRTELE DLVKVLYLER		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com