

Product Data Sheet

Recombinant Human CLSTN2 Protein

Catalog #	Source	Reactivity	Applications
CST1974	E. coli	H	
Description	Recombinant Human CLSTN2 Protein is produced by E. coli expression system and the target gene encoding His861-Asp950 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CLSTN2		
Alternative Names	CS2; Calsyntenin-2; Alcadin-gamma; Alc-gamma		
Entrez Gene	64084		
SwissProt	Q9H4D0		
Purity			
Chemical Structure	HFIQETAAK ESEMDWDDSA LTITVNPMEK HEGPGHGEDE TEGEEEEAE EEMSSSSGSD DSEEEEEEG MGRGRHGQNG ARQAQLEWDD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com