

Product Data Sheet

Recombinant Human FCRL6 Protein

Catalog #	Source	Reactivity	Applications
CST1976	E. coli	H	
Description	Recombinant Human FCRL6 Protein is produced by E. coli expression system and the target gene encoding Glu351-Arg420 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FCRL6		
Alternative Names	FCRH6; Fc receptor-like protein 6; FcR-like protein 6; FcRL6; Fc receptor homolog 6; FcRH6; IFGP6		
Entrez Gene	343413		
SwissProt	Q6DN72		
Purity			
Chemical Structure	EQCPLYANVH HQKGKDEGVV YSVVHRTSKR SEARSAEFTV GRKDSSIICA EVRCLQPSEV SSTEVMNMRSR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com