

Product Data Sheet

Recombinant Human SCGB3A1 Protein

Catalog #	Source	Reactivity	Applications
CST1982	HEK293	H	
Description	Recombinant Human SCGB3A1 Protein is produced by HEK293 expression system and the target gene encoding Phe21-Gly104 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SCGB3A1		
Alternative Names	HIN1; PNSP2; UGRP2; Secretoglobin family 3A member 1; Cytokine HIN-1; High in normal 1; Pneumo secretory protein 2; PnSP-2; Uteroglobin-related protein 2		
Entrez Gene	92304		
SwissProt	Q96QR1		
Purity			
Chemical Structure	FLVGSAPVA QPVAALESAA EAGAGTLANP LGTLNPLKLL LSSLGIPVNH LIEGSQKQVA ELGPQAVGAV KALKALLGAL TVFG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com