

Product Data Sheet

Recombinant Human TREML1 Protein

Catalog #	Source	Reactivity	Applications
CST1989	E. coli	H	
Description	Recombinant Human TREML1 Protein is produced by E. coli expression system and the target gene encoding Ser21-Ser160 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TREML1		
Alternative Names	TLT1; Trem-like transcript 1 protein; TLT-1; Triggering receptor expressed on myeloid cells-like protein 1		
Entrez Gene	340205		
SwissProt	Q86YW5		
Purity			
Chemical Structure	SLPEVLQAPV GSSILVQCHY RLQDVKAQKV WCRFLPEGCQ PLVSSAVDRR APAGRRTFLT DLGGGLLQVE MVTLQEEDAG EYGCMVDGAR GPQILHRVSL NILPPEEEEE THKIGSLAEN AFSDPAGSAN PLEPSQDEKS		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com