

Recombinant Human C1QTNF9 Protein

Catalog #	Source	Reactivity	Applications
CST1996	E. coli	H	
Description	Recombinant Human C1QTNF9 Protein is produced by E. coli expression system and the target gene encoding Pro31-Phe330 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	C1QTNF9		
Alternative Names	C1QTNF9A; Complement C1q and tumor necrosis factor-related protein 9A; Complement C1q and tumor necrosis factor-related protein 9		
Entrez Gene	338872		
SwissProt	POC862		
Purity			
Chemical Structure	PGNPGHNGLP GRDGRDGAKG DKGDAGEPGR PGSPGKDGTS GEKGERGADG KVEAKGIKGD QGSRGSPGKH GPKGLAGPMG EKGLRGETGP QGQKGNKGDV GPTGPEGPRG NIGPLGPTGL PGPMGPIGKP GPKGEAGPTG PQGEPGVRGI RGWKGDRGEK GKIGETLVLP KSAFTVGLTV LSKFPSSDMP IKFDKILYNE FNHYDTAAGK FTCHIAGVYY FTYHITVFSR NV		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com