

Product Data Sheet

Recombinant Human LRIG3 Protein

Catalog #	Source	Reactivity	Applications
CST2008	HEK293	H	
Description	Recombinant Human LRIG3 Protein is produced by HEK293 expression system and the target gene encoding Tyr871-Pro1110 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	LRIG3		
Alternative Names	LIG3; Leucine-rich repeats and immunoglobulin-like domains protein 3; LIG-3		
Entrez Gene	121227		
SwissProt	Q6UXM1		
Purity			
Chemical Structure	YVSESGSHH QFVTSSGAGF FLPQHDSSGT CHIDNSSEAD VEAATDLFLC PFLGSTGPMY LKGNVYGSDP FETYHTGCSP DPRTVLMDDHY EPSYIKKKEC YPCSHPSEES CERSFSNISW PSHVRKLLNT SYSHNEGPGM KNLCLNKSSL DFSANPEPAS VASSNSFMGT FGKALRRPHL DAYSSFGQPS DCQPRAFYLK AHSSPDLD SG SEEDGKERTD FQENHICTF KQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com