

Product Data Sheet

Recombinant Human FCRLB Protein

Catalog #	Source	Reactivity	Applications
CST2010	E. coli	H	
Description	Recombinant Human FCRLB Protein is produced by E. coli expression system and the target gene encoding Pro291-Asp370 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FCRLB		
Alternative Names	FCRL2; FCRLM2; FCRY; FREB2; Fc receptor-like B; Fc receptor homolog expressed in B-cells protein 2; FREB-2; Fc receptor-like and mucin-like protein 2; Fc receptor-like protein 2; Fc receptor-related protein Y; FcRY		
Entrez Gene	127943		
SwissProt	Q6BAA4		
Purity			
Chemical Structure	PLDPASTTAP APWAAALAPG NRPLSFRKPP VSRVPLVTS VRNTTSTGLQ FPASGAPTAG PPACAPPTPL EQSAGALKPD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com