

Product Data Sheet

Recombinant Human GSAP Protein

Catalog #	Source	Reactivity	Applications
CST2020	E. coli	H	
Description	Recombinant Human GSAP Protein is produced by E. coli expression system and the target gene encoding Leu11-Asn80 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	GSAP		
Alternative Names	PION; Gamma-secretase-activating protein; GSAP; Protein pigeon homolog		
Entrez Gene	54103		
SwissProt	A4D1B5		
Purity			
Chemical Structure	LGKDVLPWLR AQRVSEASG AGSGGADVLE NDYESLHVLN VERNGNIIYT YKDDKGNVVF GLYDCQTRQN		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com