

Product Data Sheet

Recombinant Human CDHR5 Protein

Catalog #	Source	Reactivity	Applications
CST2033	E. coli	H	
Description	Recombinant Human CDHR5 Protein is produced by E. coli expression system and the target gene encoding Ser31-Ala430 is expressed with N-terminus MBP & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CDHR5		
Alternative Names	MUCDHL; MUPCDH; Cadherin-related family member 5; Mu-protocadherin; Mucin and cadherin-like protein		
Entrez Gene	53841		
SwissProt	Q9HBB8		
Purity			
Chemical Structure	SVNKDIFEVE ENTNVTEPLV DIHVPEGQEV TLGALSTPFA FRIQGNQLFL NVTPDYEEKS LLEAQLLCQS GGTLVTLQRLV FVSVLDVNDN APEFPFKTKE IRVEEDTKVN STVIPETQLQ AEDRDKDDIL FYTLQEMTAG ASDYFSLVSV NRPALRLDRP LDFYERPNMT FWLLVRDTPG ENVEPSHTAT ATLVLVNVVPA DLRPPWFLPC TFS DG YVCIQ AQYHGAVPTG HI		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com