

Product Data Sheet

Recombinant Human ADGRG2 Protein

Catalog #	Source	Reactivity	Applications
CST2037	E. coli	H	
Description	Recombinant Human ADGRG2 Protein is produced by E. coli expression system and the target gene encoding Asp41-Pro250 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	GPR64		
Alternative Names	HE6; TM7LN2; G-protein coupled receptor 64; Human epididymis-specific protein 6; He6		
Entrez Gene	10149		
SwissProt	Q8IZP9		
Purity			
Chemical Structure	DTDNSSLSP PAKLSVVSFA PSSNGTPEVE TTSLNDVTLS LLPSNETEKT KITIVKTFNA SGVKPQRNIC NLSSICNDSA FFRGEIMFYQ DKESTVPQNN HITNGTLTGVL LSLSELKRSE LNKTLQTLSE TYFIMCATAE AQSTLNCTFT IKLNNTMNAC AVIAALERVK IRPMEHCCCS VRIPCPSSPE ELEKLQCDLQ DPIVCLADHP		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com