

# Product Data Sheet

## Recombinant Human FST Protein

| Catalog #                 | Source                                                                                                                                                                                                                                                                | Reactivity | Applications |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|--------------|
| CST2049                   | HEK293                                                                                                                                                                                                                                                                | H          |              |
| <b>Description</b>        | Recombinant Human FST Protein is produced by HEK293 expression system and the target gene encoding Asn31-Ser320 is expressed with C-terminus His tag.                                                                                                                 |            |              |
| <b>Form</b>               | Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.                                                                                                                                                                                                                  |            |              |
| <b>Gene Symbol</b>        | FST                                                                                                                                                                                                                                                                   |            |              |
| <b>Alternative Names</b>  | Follistatin; FS; Activin-binding protein                                                                                                                                                                                                                              |            |              |
| <b>Entrez Gene</b>        | 10468                                                                                                                                                                                                                                                                 |            |              |
| <b>SwissProt</b>          | P19883                                                                                                                                                                                                                                                                |            |              |
| <b>Purity</b>             |                                                                                                                                                                                                                                                                       |            |              |
| <b>Chemical Structure</b> | NCWLRQAKNG RCQVLYKTEL SKEECSTGR LTSWTEEDV NDNTLFKWMI FNGGAPNCIP<br>CKETCENVDC GPGKKCRMNK KMKPRVCAP DCSNITWKGP VCGLDGKTYR NECALLKARC<br>KEQPELEVQY QGRCKKTCRD VFCPGSSTCV VDQTNNAVCV TCNRICPEPA SSEQYLCGND<br>GVTYSSACHL RKATCLLGRS IGLAYEGKCI KAKSCEDIQC TGGKKCLWDF KV |            |              |
| <b>Quality Control</b>    |                                                                                                                                                                                                                                                                       |            |              |
| <b>Directions for Use</b> | Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.                                                                                                                                                     |            |              |
| <b>Storage/Stability</b>  | Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.                                                                                                                      |            |              |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com