

Product Data Sheet

Recombinant Human FCRL1 Protein

Catalog #	Source	Reactivity	Applications
CST2067	E. coli	H	
Description	Recombinant Human FCRL1 Protein is produced by E. coli expression system and the target gene encoding Pro351-Ile420 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FCRL1		
Alternative Names	FCRH1; IFGP1; IRTA5; Fc receptor-like protein 1; FcR-like protein 1; FcRL1; Fc receptor homolog 1; FcRH1; IFGP family protein 1; hIFGP1; Immune receptor translocation-associated protein 5; CD307a		
Entrez Gene	115350		
SwissProt	Q96LA6		
Purity			
Chemical Structure	PQEFTYLNSP TPGQLQPIYE NVNVVSGDEV YSLAYYNQPE QESVAAETLG THMEDKVSLE IYSRLRKANI		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com