

Product Data Sheet

Recombinant Human DCBLD2 Protein

Catalog #	Source	Reactivity	Applications
CST2077	E. coli	H	
Description	Recombinant Human DCBLD2 Protein is produced by E. coli expression system and the target gene encoding Gly321-Thr520 is expressed with N-terminus SUMO & N-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	DCBLD2		
Alternative Names	CLCP1; ESDN; Discoidin, CUB and LCCL domain-containing protein 2; CUB, LCCL and coagulation factor V/VIII-homology domains protein 1; Endothelial and smooth muscle cell-derived neuropilin-like protein		
Entrez Gene	131566		
SwissProt	Q96PD2		
Purity			
Chemical Structure	GQENSWKPKK ARLKKPGPPW AAFATDEYQW LQIDLNKEKK ITGIITTGST MVEHNYVSA YRILYSDDGQ KWTVYREPGV EQDKIFQGNK DYHQDVRNNF LPPIARFIR VNPTQWQQKI AMKMELLGCQ FIPKGRPPKL TQPPPPRNSN DLKNTTAPPK IAKGRAPKFT QPLQPRSSNE FPAQTEQTTA SPDIRNTTVT		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com