

Product Data Sheet

Recombinant Human FOXJ3 Protein

Catalog #	Source	Reactivity	Applications
CST2105	E. coli	H	
Description	Recombinant Human FOXJ3 Protein is produced by E. coli expression system and the target gene encoding Ser221-Ala470 is expressed with N-terminus MBP & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FOXJ3		
Alternative Names	KIAA1041; Forkhead box protein J3		
Entrez Gene	22887		
SwissProt	Q9UPW0		
Purity			
Chemical Structure	SDSPRSSLNN SLSDQSLASV NLNSVGSVHS YTPVTSHPEVSQSLTPQQQ PQYNLPERDK QLLFSEYNFE DLSASFRSLY KSVFEQSLSQ QGLMNIPSES SQQSHTSCTY QHSPSSTVST HPHSNQSSLS NSHGSGLNNTT GSNSVAQVSL SHPQMHTQPS PHPPHRPHGL PQHPQRSPHP AHPQQHSQQL QSPHPQHPSP HQHIQHHPNH QHQTLT HQAP PPPQQVSCNS GV		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com