

Recombinant Human FCRL3 Protein

Catalog #	Source	Reactivity	Applications
CST2156	E. coli	H	
Description	Recombinant Human FCRL3 Protein is produced by E. coli expression system and the target gene encoding Asp631-Leu730 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FCRL3		
Alternative Names	FCRH3; IFGP3; IRTA3; SPAP2; Fc receptor-like protein 3; FcR-like protein 3; FcRL3; Fc receptor homolog 3; FcRH3; IFGP family protein 3; hIFGP3; Immune receptor translocation-associated protein 3; SH2 domain-containing phosphatase anchor protein 2; CD307c		
Entrez Gene	115352		
SwissProt	Q96P31		
Purity			
Chemical Structure	DPQEPTHSKP LAPMELEPMY SNVNP GDSNP IYSQIWSIQH TKENSANCPM MHQEHEELTV LYSELKKTHP DDSAGEASSR GRAHEEDDEE NYENVPRVLL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com