

Product Data Sheet

Recombinant Human PIBF1 Protein

Catalog #	Source	Reactivity	Applications
CST2163	E. coli	H	
Description	Recombinant Human PIBF1 Protein is produced by E. coli expression system and the target gene encoding Lys11-Glu550 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	PIBF1		
Alternative Names	C13orf24; PIBF; Progesterone-induced-blocking factor 1		
Entrez Gene	10464		
SwissProt	Q8WXW3		
Purity			
Chemical Structure	KVNISSSLES EDISLETTVP TDDISSEER EGKVRITRQL IERKELLHNI QLLKIELSQK TMMIDNLKVD YLTKEELEE KLNDALHQKQ LLTLRLDNQL AFQQKDASKY QELMKQEMET ILLRQKQLEE TNLQLREKAG DVRRNLRDFE LTEEQYIKLK AFPEDQLSIP EYVSVRFYEL VNPLRKEICE LQVKKNILAE ELSTNKNQLK QLTETYEEDR KNYSEVQIRC QR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com