

Product Data Sheet

Recombinant Human KIRREL3 Protein

Catalog #	Source	Reactivity	Applications
CST2169	E. coli	H	
Description	Recombinant Human KIRREL3 Protein is produced by E. coli expression system and the target gene encoding Val571-Ser680 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	KIRREL3		
Alternative Names	KIAA1867; NEPH2; Kin of IRRE-like protein 3; Kin of irregular chiasm-like protein 3; Nephrin-like protein 2		
Entrez Gene	84623		
SwissProt	Q8IZU9		
Purity			
Chemical Structure	VSAKNDIRVE IVHKEPASGR EGEEHSTIKQ LMMDRGEFQQ DSVLKQLEVL KEEKEFQNL KDPTNGYYSV NTFKEHHSTP TISLSSCQPD LRPAGKQRVP TGMSFTNIYS		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com