

Product Data Sheet

Recombinant Human MCEMP1 Protein

Catalog #	Source	Reactivity	Applications
CST2171	E. coli	H	
Description	Recombinant Human MCEMP1 Protein is produced by E. coli expression system and the target gene encoding Glu11-Tyr80 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	MCEMP1		
Alternative Names	C19orf59; Mast cell-expressed membrane protein 1		
Entrez Gene	199675		
SwissProt	Q8IX19		
Purity			
Chemical Structure	EVKMQAPAFR DKKQGVSAKN QGAHDPDYEN ITLAFKNQDH AKGGHSRPTS QVPAQCRPPS DSTQVPCWLY		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com