

Product Data Sheet

Recombinant Human APBB1IP Protein

Catalog #	Source	Reactivity	Applications
CST2173	E. coli	H	
Description	Recombinant Human APBB1IP Protein is produced by E. coli expression system and the target gene encoding Arg81-Lys260 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	APBB1IP		
Alternative Names	PREL1; RARP1; RIAM; Amyloid beta A4 precursor protein-binding family B member 1-interacting protein; APBB1-interacting protein 1; Proline-rich EVH1 ligand 1; PREL-1; Proline-rich protein 73; Rap1-GTP-interacting adapter molecule; RIAM; Retinoic acid-respo		
Entrez Gene	54518		
SwissProt	Q7Z5R6		
Purity			
Chemical Structure	RTIQAQKESL QNQHHSASLQ ASIFSGAASL GYGTVNAATG ISQYEDDLPP PPADPVLDLP LPPPPPEPLS QEEEEAQAKA DKIKLALAKL KEAKVKKLTV KVHMNDNSTK SLMVDERQLA RDVLDNLFK THCDCNVDCW LYEIPELQI ERFFEDHENV VEVLSDWTRD TENKILFLEK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com