

Product Data Sheet

Recombinant Human DAB2 Protein

Catalog #	Source	Reactivity	Applications
CST2199	E. coli	H	
Description	Recombinant Human DAB2 Protein is produced by E. coli expression system and the target gene encoding Asn11-Ile260 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	DAB2		
Alternative Names	DOC2; Disabled homolog 2; DOC-2; Differentially-expressed protein 2		
Entrez Gene	1601		
SwissProt	P98082		
Purity			
Chemical Structure	NGQPDQQAAP KAPSKKEKKK GPEKTDEYLL ARFKGDGVKY KAKLIGIDDV PDARGDKMSQ DSMMKLKGMA AAGRSQGQHK QRIWVNISLS GIKIIDEKTG VIEHEHPV NK ISFIARDVTD NRAFGYVCGG EGQHQQFAIK TGQQAELVV DLKDLFQVIY NVKKKEEEEK KIEEASKAVE NGSEALMILD DQTNKLKSGV DQMDLFGDMS TPPDLNSPTE SKDILLVDLN SE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com