

Product Data Sheet

Recombinant Human RAP1A Protein

Catalog #	Source	Reactivity	Applications
CST2201	E. coli	H	
Description	Recombinant Human RAP1A Protein is produced by E. coli expression system and the target gene encoding Lys31-Gln130 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	RAP1A		
Alternative Names	KREV1; Ras-related protein Rap-1A; C21KG; G-22K; GTP-binding protein smg p21A; Ras-related protein Krev-1		
Entrez Gene	5906		
SwissProt	P62834		
Purity			
Chemical Structure	KYDPTIEDSY RKQVEVDCQQ CMLEILD TAG TEQFTAMRDL YMKNGQG FAL VYSIT AQSTF NDLQDLREQI LRVKDTEDVP MILVG NKCDL EDERVVGKEQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com