

Product Data Sheet

Recombinant Human TXLNA Protein

Catalog #	Source	Reactivity	Applications
CST2209	E. coli	H	
Description	Recombinant Human TXLNA Protein is produced by E. coli expression system and the target gene encoding Ala11-Gln160 is expressed with N-terminus MBP & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TXLNA		
Alternative Names	TXLN; Alpha-taxilin		
Entrez Gene	200081		
SwissProt	P40222		
Purity			
Chemical Structure	AKQSNPKSSP GQPEAGPEGA QERPSQAAPA VEAEGPGSSQ APRKPEGAQA RTAQSGALRD VSEELSRQLE DILSTYCVDN NQGGPGEDGA QGEPAPEDA EKSRTYVARN GEPEPTPVVN GEKEPSKGDP NTEEIRQSDE VGDRDHRRPQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com