

Product Data Sheet

Recombinant Human GLI2 Protein

Catalog #	Source	Reactivity	Applications
CST2225	E. coli	H	
Description	Recombinant Human GLI2 Protein is produced by E. coli expression system and the target gene encoding Glu11-Ser90 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	GLI2		
Alternative Names	THP; Zinc finger protein GLI2; Tax helper protein		
Entrez Gene	2736		
SwissProt	P10070		
Purity			
Chemical Structure	EKQEAKSGIL EAAGFPDPGK KASPLVAAA AAAAVAAQGV PQHLLPPFHA PLPIDMRHQE GRYHYEPHSV HGVHGPPALS		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com