

Product Data Sheet

Recombinant Human GALNT10 Protein

Catalog #	Source	Reactivity	Applications
CST2250	E. coli	H	
Description	Recombinant Human GALNT10 Protein is produced by E. coli expression system and the target gene encoding Gly41-Glu110 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	GALNT10		
Alternative Names	Polypeptide N-acetylgalactosaminyltransferase 10; Polypeptide GalNAc transferase 10; GalNAc-T10; pp-GaNTase 10; Protein-UDP acetylgalactosaminyltransferase 10; UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 10		
Entrez Gene	55568		
SwissProt	Q86SR1		
Purity			
Chemical Structure	GGSGAAVAPA AGQGSHSRQK KTFFLGDGQK LKDWHDKEAI RRDAQRVGNG EQGRPYPM TD AERVDQAYRE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com