

Product Data Sheet

Recombinant Human MAD1L1 Protein

Catalog #	Source	Reactivity	Applications
CST2251	E. coli	H	
Description	Recombinant Human MAD1L1 Protein is produced by E. coli expression system and the target gene encoding Leu11-Asp350 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	MAD1L1		
Alternative Names	MAD1; TXBP181; Mitotic spindle assembly checkpoint protein MAD1; Mitotic arrest deficient 1-like protein 1; MAD1-like protein 1; Mitotic checkpoint MAD1 protein homolog; HsMAD1; hMAD1; Tax-binding protein 181		
Entrez Gene	8379		
SwissProt	Q9Y6D9		
Purity			
Chemical Structure	LSTLRSLNNF ISQRVEGGSG LDISTSAPGS LQMQYQQSMQ LEERAEQIRS KSHLIQVERE KMQMELSHKR ARVELERAAS TSARNYEREV DRNQELLTRI RQLQEREAGA EEKMQEQLER NRQCQQNLDA ASKRLREKED SLAQAGETIN ALKGRISELQ WSVMDQEMRV KRLESEKQEL QEQLDLQHKK CQEANQKIQE LQASQEARAD HEQQIKDLEQ KLSLQEQDAA IV		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com