

Product Data Sheet

Recombinant Human CDON Protein

Catalog #	Source	Reactivity	Applications
CST2253	E. coli	H	
Description	Recombinant Human CDON Protein is produced by E. coli expression system and the target gene encoding Phe31-Arg960 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CDON		
Alternative Names	CDO; Cell adhesion molecule-related/down-regulated by oncogenes		
Entrez Gene	50937		
SwissProt	Q4KMG0		
Purity			
Chemical Structure	FTSEPLSAVQ KLGGPVVLHC SAQPVTTTRIS WLHNGKTLDG NLEHVKIHQG TLTILSLNSS LLGYQCLAN NSIGAIVSGP ATVSVAVLGD FGSSTKHVIT AEEKSAGFIG CRVPESNPKA EVRYKIRGKW LEHSTENYLI LPSGNLQILN VSLEDKGSYK CAAYNPVTHQ LKVEPIGRKL LVSRPSSDDV HILHPTHSAQ LAVLSRSPVT LECVVSGVPA PQVYWLKDGQ DI		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com