

Product Data Sheet

Recombinant Human GFRA2 Protein

Catalog #	Source	Reactivity	Applications
CST2256	E. coli	H	
Description	Recombinant Human GFRA2 Protein is produced by E. coli expression system and the target gene encoding Asp341-Leu410 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	GFRA2		
Alternative Names	GDNFRB; RETL2; TRNR2; GDNF family receptor alpha-2; GDNF receptor alpha-2; GDNFR-alpha-2; GFR-alpha-2; GDNF receptor beta; GDNFR-beta; Neurturin receptor alpha; NRTNR-alpha; NTN-2; RET ligand 2; TGF-beta-related neurotrophic factor receptor 2		
Entrez Gene	2675		
SwissProt	O00451		
Purity			
Chemical Structure	DFTENPCLRN AIQAFGNGTD VNVSPKGPSF QATQAPRVEK TPSLPDDLSD STSLGTSVIT TCTSVQEQGL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com