

Product Data Sheet

Recombinant Human TRAF3 Protein

Catalog #	Source	Reactivity	Applications
CST2262	E. coli	H	
Description	Recombinant Human TRAF3 Protein is produced by E. coli expression system and the target gene encoding His251-Ala400 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TRAF3		
Alternative Names	CAP1; CRAF1; TNF receptor-associated factor 3; CAP-1; CD40 receptor-associated factor 1; CRAF1; CD40-binding protein; CD40BP; LMP1-associated protein 1; LAP1		
Entrez Gene	7187		
SwissProt	Q13114		
Purity			
Chemical Structure	HEASSAVQHV NLLKEWSNSL EKKVSLQNE SVEKNKSIQS LHNQICSFEI EIERQKEMLR NNESKILHLQ RVIDSQAELK KELDKEIRPF RQNWEEADSM KSSVESLQNR VTELESVDKS AGQVARNTGL LESQLSRHDQ MLSVHDIRLA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com