

Product Data Sheet

Recombinant Human SIGLEC9 Protein

Catalog #	Source	Reactivity	Applications
CST2283	E. coli	H	
Description	Recombinant Human SIGLEC9 Protein is produced by E. coli expression system and the target gene encoding Ala381-Lys460 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SIGLEC9		
Alternative Names	Sialic acid-binding Ig-like lectin 9; Siglec-9; CDw329; Protein FOAP-9; CD329		
Entrez Gene	27180		
SwissProt	Q9Y336		
Purity			
Chemical Structure	AAGVGDTGIE DANA VRGSAS QGPLTEPWAE DSPPDQPPPA SARSSVGEGE LQYASLSFQM VKPWDSRGQE ATDTEYSEIK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com