

Product Data Sheet

Recombinant Human SENP2 Protein

Catalog #	Source	Reactivity	Applications
CST2288	E. coli	H	
Description	Recombinant Human SENP2 Protein is produced by E. coli expression system and the target gene encoding Asn91-Ser390 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SEN2		
Alternative Names	KIAA1331; Sentrin-specific protease 2; Axam2; SMT3-specific isopeptidase 2; Smt3ip2; Sentrin/SUMO-specific protease SENP2		
Entrez Gene	59343		
SwissProt	Q9HC62		
Purity			
Chemical Structure	NVAPSGEVFS NSSCELTS GSWNNMLKL NKSPNGISDY PKIRVTVTRD QPRRVLPSTG FTLNSEGCNR RPPGRRHSG NPESLMWKP QEQAVTEMIS EESGKGLRRP HCTVEEGVQK EEREKYRKLL ERLKESGHGN SVCPVTSNYH SSQRSQMDTL KTKGWGEEQN HGVKTTQFVP KQYRLVETRG PLCSLRSEKR CSKGKITDTE TMVGIRFENE SRRGYQLEPD LS		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com