

Product Data Sheet

Recombinant Human Extracellular Matrix Protein 1 Protein

Catalog #	Source	Reactivity	Applications
CST2295	E. coli	H	
Description	Recombinant Human Extracellular Matrix Protein 1 Protein is produced by E. coli expression system and the target gene encoding Leu51-His150 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ECM1		
Alternative Names	Extracellular matrix protein 1; Secretory component p85		
Entrez Gene	1893		
SwissProt	Q16610		
Purity			
Chemical Structure	LSRSLPMDHP DSSQHGPPE GQSQVQPPPS QEATPLQKEK LLPAQLPAEK EVGPPLPQEA VPLQKELPSL QHPNEQKEGT PAPFGDQSHP EPESWNAAQH		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com