

Product Data Sheet

Recombinant Human NOMO1 Protein

Catalog #	Source	Reactivity	Applications
CST2298	E. coli	H	
Description	Recombinant Human NOMO1 Protein is produced by E. coli expression system and the target gene encoding Ser701-Phe880 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	NOMO1		
Alternative Names	PM5; Nodal modulator 1; pM5 protein		
Entrez Gene	23420		
SwissProt	Q15155		
Purity			
Chemical Structure	SVQELRREQQ LAEIEARRQE REKNGNEEGE ERMTKPPVQE MVDELQGPFS YDFSYWARGSG EKITVTPSSK ELLFYPPSME AVVSGESCPG KLIEIHGKAG LFLEGQIHPE LEGVEIVISE KGASSPLITV FTDDKGAYSV GPLHSDLEYT VTSQKEGYVL TAVEGTIGDF KAYALAGVSF		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com