

Product Data Sheet

Recombinant Human ZBTB17 Protein

Catalog #	Source	Reactivity	Applications
CST2302	E. coli	H	
Description	Recombinant Human ZBTB17 Protein is produced by E. coli expression system and the target gene encoding Gly651-Asp750 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ZBTB17		
Alternative Names	MIZ1; ZNF151; ZNF60; Zinc finger and BTB domain-containing protein 17; Myc-interacting zinc finger protein 1; Miz-1; Zinc finger protein 151; Zinc finger protein 60		
Entrez Gene	7709		
SwissProt	Q13105		
Purity			
Chemical Structure	GSEVSVVTVD DMVTLATEAL AATAVTQLTV VPVGAAVTAD ETEVLKAEIS KAVKQVQEED PNTHILYACD SCGDKFLDAN SLAQHVRIHT AQALVMFQTD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com