

## Recombinant Human STK4 Protein

| Catalog #                 | Source                                                                                                                                                                                                        | Reactivity | Applications |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|--------------|
| CST2303                   | E. coli                                                                                                                                                                                                       | H          |              |
| <b>Description</b>        | Recombinant Human STK4 Protein is produced by E. coli expression system and the target gene encoding Lys301-Lys480 is expressed with C-terminus His tag.                                                      |            |              |
| <b>Form</b>               | Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.                                                                                                                                                          |            |              |
| <b>Gene Symbol</b>        | STK4                                                                                                                                                                                                          |            |              |
| <b>Alternative Names</b>  | KRS2; MST1; Serine/threonine-protein kinase 4; Mammalian STE20-like protein kinase 1; MST-1; STE20-like kinase MST1; Serine/threonine-protein kinase Krs-2                                                    |            |              |
| <b>Entrez Gene</b>        | 6789                                                                                                                                                                                                          |            |              |
| <b>SwissProt</b>          | Q13043                                                                                                                                                                                                        |            |              |
| <b>Purity</b>             |                                                                                                                                                                                                               |            |              |
| <b>Chemical Structure</b> | KLKRQESQQR EVDQDDEENS EEDEMDSGTM VRAVGDEMGT VRVASTMTDG<br>ANTMIEHDDT LPSQLGTMVI NAEDEEEEGT MKRRDETMQP AKPSFLEYFE QKEKENQINS<br>FGKSVPGPLK NSSDWKIPQD GDYEFKSWT VEDLQKRLLA LDPMMEQEIE EIRQKYQSKR<br>QPILDAIEAK |            |              |
| <b>Quality Control</b>    |                                                                                                                                                                                                               |            |              |
| <b>Directions for Use</b> | Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.                                                                                             |            |              |
| <b>Storage/Stability</b>  | Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.                                                              |            |              |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)