

Product Data Sheet

Recombinant Human ANXA11 Protein

Catalog #	Source	Reactivity	Applications
CST2314	E. coli	H	
Description	Recombinant Human ANXA11 Protein is produced by E. coli expression system and the target gene encoding Asp291-Asn370 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ANXA11		
Alternative Names	ANX11; Annexin A11; 56 kDa autoantigen; Annexin XI; Annexin-11; Calcyclin-associated annexin 50; CAP-50		
Entrez Gene	311		
SwissProt	P50995		
Purity			
Chemical Structure	DEACLEILA SRSNEHIREL NRAYKAEFFK TLEEAI RSDT SGHFQRL LIS LSQGNRDEST NVDM SLAQ RD AQELYAAGEN		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com