

Product Data Sheet

Recombinant Human MST1 Protein

Catalog #	Source	Reactivity	Applications
CST2326	E. coli	H	
Description	Recombinant Human MST1 Protein is produced by E. coli expression system and the target gene encoding Tyr371-Lys470 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	MST1		
Alternative Names	D3F15S2; DNF15S2; HGFL; Hepatocyte growth factor-like protein; Macrophage stimulatory protein; Macrophage-stimulating protein; MSP		
Entrez Gene	4485		
SwissProt	P26927		
Purity			
Chemical Structure	YHGAGEQYRG TVSKTRKGVQ CQRWSAETPH KPQFTFTSEP HAQLEENFCR NPDGDSHGPDW CYTMDPRTPF DYCALRRCAD DQPPSILDPP DQVQFEKCGK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com